E

1- DESCRIPCIÓN

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


EN

1- DESCRIPTION

This brace is designed to stabilize the motion of the knee joint, keeping it aligned during flexion and extension movements. Constructed entirely from breathable materials, the next-to-skin layer is made from anti-bacterial and anti-allergy fabric. Its complex structure, equipped with patella stabilizer, offers comfortable containment of the patellar compartment, useful in cases of associated conditions (medial or lateral ligaments, anterior or posterior cruciate ligament and patellar instability). Two patented, multicentric, adjustable aluminum rods are included to control flexion and extension quickly and easily. The rods are equipped with bearings supporting the condyles and rigid tibial and femoral supports to control rotation. This brace conforms to the requirements of policy 93/42/EEC re. medical devices and Italian legislation (DL 24 February 1997, no. 46).

2- INDICATIONS

A knee orthosis with adjustable knee joints used for:

- Patellar chondropathy (osteochondral repairs)
- Retropatellar osteoarthritis
- Post-operative and post-traumatic irritations including rehabilitation of: MCL, LCL, PCL, ACL (Sprains, Meniscus injuries, Lateral, Anterior, Medial & Posterior tilt, etc.)
- Post lateral release surgery to limit the extension and flexion of the joint
- Patella and knee instabilities (Patella Dislocation, Patella Subluxation, Patella Alta, Patella stabilization, Patella tendon repairs)
- General knee pain
- Mild Chondromalacia
- Condylar Fractures
- Tibial Plateau Fractures

2-1 INSTRUCTIONS

Apply brace directly onto skin following these instructions:

1. Open locking system Velcro;
2. Pull on the brace (keeping the label on the top). Position in such a way that the (C) pressure pads surround the patella;
3. Close the three fasteners by pulling on both parts of the Velcro of each fastener at the same time, in the following order; (A) Middle, (B) Bottom, (C) Top.
4. Apply the Velcro panel to the back of the thigh. This application will conform to the shape of the thigh.
5. After applying the Velcro panel, tighten the straps "C" using the normal locking technique one time more.
6. Metal Rod (ROM - Range Of Motion). For Rod extension and flexion adjustment, the measurements must be indicated by the attending physician. For adjustment, operate as follows:
Open the side zippers to reveal the rigid internal supports; two Colored tabs appear, the red for "Flexion" and grey for "Extension". Pull the Flexion tab so that the loop wraps around the hook showing the number of desired degree. Do the same with the Extension tab. Repeat the same operation on the other rod then make sure that both rods have the same adjustments.

3- SIDE EFFECTS AND CONTRA-INDICATIONS

Assuming that the products are used according to instructions and given that they are non invasive, there are no side effects affecting the organism as a whole. All braces and support, if overly tight, can cause discomfort due to compression. If you observe any change in your general state. During a prolonged period of inactivity, it is recommended to reduce the tension on your device until you resume normal activities where you will reset the tension to the original level recommended to you.
Any modifications of these products will void the warranty along with the denial of any responsibilities on our behalf.

4- MAINTENANCE AND CARE

- No special maintenance is required.
- Attach the Velcro pads firmly to the surface.
- Wash gently at lukewarm temperature (machine wash or hand wash).
- Do not use any substances containing chlorine.
- Dry at low setting.
- Washing does not damage the product and it actually restores the elasticity.

This product is not biodegradable.
AT THE END OF ITS USABLE LIFE, DISCARD IT ALONG WITH YOUR NORMAL WASTE.

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